

**Fast Pap
05 - 12055**

IVD In-vitro diagnostic medical device **CE**
CND Code: W01030708

number of tests 300-500
procedure time 3 minutes
product validity 1 year
storage temperature 15-25°C

Expected aim	Product for the preparation of cyto-histological samples for optical microscopy.																						
Application	Gynecological exfoliative cytology, body fluids, fine needle aspiration cytology, intraoperative cytology. Liquid based cytology.																						
Principle	A highly selective blue nuclear stain, Harris' hematoxylin, is combined with IW3 polychromic mixture, a subtle cytoplasmic stain which differentiates cyanophil cells from eosinophil ones.																						
Method	<table border="0"> <tr><td>1) Distilled water</td><td>10 seconds</td></tr> <tr><td>2) Nuclear stain - Gill 3 Hematoxylin</td><td>45 seconds</td></tr> <tr><td>3) Colour change in tap water</td><td>15 seconds</td></tr> <tr><td>4) Dehydrating agent - Unyhol</td><td>5 seconds</td></tr> <tr><td>5) Cytoplasmatic stain - IW3 staining solution</td><td>1 minute</td></tr> <tr><td>6) Dehydrating agent - Unyhol</td><td>5 seconds</td></tr> <tr><td>7) Dehydrating agent - Unyhol</td><td>5 seconds</td></tr> <tr><td>8) Dehydrating agent - Unyhol</td><td>5 seconds</td></tr> <tr><td>9) Clearing agent - Bio-Clear</td><td>5 seconds</td></tr> <tr><td>10) Clearing agent - Bio-Clear</td><td>5 seconds</td></tr> <tr><td>11) Mount in xylene based mounting medium</td><td></td></tr> </table>	1) Distilled water	10 seconds	2) Nuclear stain - Gill 3 Hematoxylin	45 seconds	3) Colour change in tap water	15 seconds	4) Dehydrating agent - Unyhol	5 seconds	5) Cytoplasmatic stain - IW3 staining solution	1 minute	6) Dehydrating agent - Unyhol	5 seconds	7) Dehydrating agent - Unyhol	5 seconds	8) Dehydrating agent - Unyhol	5 seconds	9) Clearing agent - Bio-Clear	5 seconds	10) Clearing agent - Bio-Clear	5 seconds	11) Mount in xylene based mounting medium	
1) Distilled water	10 seconds																						
2) Nuclear stain - Gill 3 Hematoxylin	45 seconds																						
3) Colour change in tap water	15 seconds																						
4) Dehydrating agent - Unyhol	5 seconds																						
5) Cytoplasmatic stain - IW3 staining solution	1 minute																						
6) Dehydrating agent - Unyhol	5 seconds																						
7) Dehydrating agent - Unyhol	5 seconds																						
8) Dehydrating agent - Unyhol	5 seconds																						
9) Clearing agent - Bio-Clear	5 seconds																						
10) Clearing agent - Bio-Clear	5 seconds																						
11) Mount in xylene based mounting medium																							
Results	<table border="0"> <tr><td>Nuclei.....</td><td>blue-violet</td></tr> <tr><td>Cyanophilic cytoplasm.....</td><td>blue-green</td></tr> <tr><td>Eosinophilic cytoplasm.....</td><td>pink</td></tr> <tr><td>Keratinized cytoplasm.....</td><td>orange-pink</td></tr> <tr><td>Microorganisms.....</td><td>blue-violet</td></tr> <tr><td>Erythrocytes.....</td><td>red-brown</td></tr> </table>	Nuclei.....	blue-violet	Cyanophilic cytoplasm.....	blue-green	Eosinophilic cytoplasm.....	pink	Keratinized cytoplasm.....	orange-pink	Microorganisms.....	blue-violet	Erythrocytes.....	red-brown										
Nuclei.....	blue-violet																						
Cyanophilic cytoplasm.....	blue-green																						
Eosinophilic cytoplasm.....	pink																						
Keratinized cytoplasm.....	orange-pink																						
Microorganisms.....	blue-violet																						
Erythrocytes.....	red-brown																						
Reagents	<table border="0"> <tr><td>A) Nuclear stain - Gill 3 Hematoxylin</td><td>1 x 500 ml</td></tr> <tr><td>B) Cytoplasmatic stain - IW3 staining solution</td><td>1 x 500 ml</td></tr> <tr><td>C) Dehydrating agent - Unyhol (2-propanol 75% – ethanol 25% mixture)</td><td>1 x 2500 ml</td></tr> <tr><td>D) Clearing agent - Bio-Clear</td><td>1 x 2500 ml</td></tr> </table>	A) Nuclear stain - Gill 3 Hematoxylin	1 x 500 ml	B) Cytoplasmatic stain - IW3 staining solution	1 x 500 ml	C) Dehydrating agent - Unyhol (2-propanol 75% – ethanol 25% mixture)	1 x 2500 ml	D) Clearing agent - Bio-Clear	1 x 2500 ml														
A) Nuclear stain - Gill 3 Hematoxylin	1 x 500 ml																						
B) Cytoplasmatic stain - IW3 staining solution	1 x 500 ml																						
C) Dehydrating agent - Unyhol (2-propanol 75% – ethanol 25% mixture)	1 x 2500 ml																						
D) Clearing agent - Bio-Clear	1 x 2500 ml																						
Warning and precaution	<p>The product must be used exclusively by specialized technical operators. Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental release.</p> <p>Do not use if the primary container is damaged.</p>																						
Storage	Store the preparation at room temperature. Keep the containers tightly closed.																						
Stability	After the first opening, the product is usable until the expiry date, if correctly stored.																						
Disposal	Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.																						

